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
TOTAL 5.407,00

## FISH in myeloma

EasySep<sup>®</sup> is a powerful immunomagnetic cell selection procedure that uses specific antibodies and tiny FACS-compatible magnetic nanoparticles in a column-free magnetic system. Cells are targeted for selection or depletion using monoclonal antibodies directed against specific cell surface antigens. These labeled cells are then crosslinked to EasySep<sup>®</sup> magnetic nanoparticles in a standard FACS tube. The tube is then placed directly in the specially designed EasySep<sup>®</sup> magnet. This handheld magnet is gently inverted to remove the cells that are not bound to the magnetic nanoparticles. The unique EasySep<sup>®</sup> magnet generates a high gradient magnetic field in the interior cavity that is strong enough to separate cells labeled with EasySep<sup>®</sup> nanoparticles without the additional magnetic field gradients provided by a column matrix.

For Positive Cell Selection: Cells of interest remain in the tube after pouring off the cells that are not bound to the magnetic nanoparticles. To obtain a highly pure population, these cells are rinsed twice or more, then easily collected by removing the tube from the magnet.